

Installation and Software Usage Guide

Pre-installation:

install XWin Server (for Windows) <http://x.cygwin.com/devel/server/>



Xterm for Linux and Mac OS <http://invisible-island.net/xterm/>

or any other applications that support X11 forwarding.

Connecting to Cluster:

SSH -X useraccount@eagle.bmi.stonybrook.edu

then

SSH -X nfs001

Download and Install on cluster:

1. eclipse-jee-luna-SR1-linux-gtk-x86_64.tar.gz (Eclipse JEE is required to run with Apache) from http://www.eclipse.org/downloads/download.php?file=/technology/epp/downloads/release/luna/SR1/eclipse-jee-luna-SR1-linux-gtk-x86_64.tar.gz
2. Apache 8 from <http://supergsego.com/apache/tomcat/tomcat-8/v8.0.15/bin/apache-tomcat-8.0.15.tar.gz>
3. Firefox Nightly from https://ftp.mozilla.org/pub/mozilla.org/firefox/nightly/latest-trunk/firefox-37.0a1.en-US.linux-x86_64.tar.bz2

for tar.gz use: tar xzf eclipse-jee-luna-SR1-linux-gtk-x86_64.tar.gz

and : tar xzf apache-tomcat-8.0.15.tar.gz


for tar.be2 use: tar xjf firefox-37.0a1.en-US.linux-x86_64.tar.bz2

Launch and configuration:


go to Eclipse directory then ./eclipse & to launch eclipse jee

and then go to FireFox directory ./firefox & to launch firefox


In Eclipse import the project (TptWebsite). There may be library version unmatched errors. To solve this, got to Configure Build Path, remove unmatched JRE System Library, use the default JRE library used by Eclipse.




Then run as-> Run on Server




For the first time use, set up new Tomcat V8.0 Sever



Move TptWebsite to Configured.




Also modify server.xml under Servers->Tomcat v8.0 Server at localhost-config



Change port. Default is 8080 (We use 8888) If there are any other conflicts, change ports to any unused port.

```
-->
<Connector connectionTimeout="20000" port="8888" protocol="HTTP/1.1" redirectPort="8443"/>
<!-- & "Connector" using the shared thread pool-->
```

Then start Server



In FireFox, go to http://localhost:port_number/TptWebsite/
(<http://localhost:8888/TptWebsite/> in our case)

Application

For Basic:

User can select search condition on Project, TSS, Participant, Sample and Vial

The screenshot shows a web application interface for 'Project 4'. It has a header with 'About', 'Services', and 'Contact' links. The main content area is titled 'Please select search condition'. Below this title, there are five dropdown menus for selecting search conditions:

- Project: TCGA
- TSS: 06
- Participant: 0208
- Sample: 01
- Vial: Z

User can also add features to search. Add Condition, then select specific feature and operator, type value. There are range intervals for each feature.

Project 4

AboutServicesContact

Feature

Add Condition

FSD1

0.1

FSD1 Range is [0.010303, 0.852794]

SEARCH

BACK TO HOMEPAGE

[Home](#) • [About](#) • [Services](#) • [Contact](#)

Copyright © Song 2014. All Rights Reserved

If there are more conditions. Keep adding conditions.

Feature

Add Condition

FSD1

0.5

FSD1 Range is [0.010303, 0.852794]

FSD2

0.2

FSD2 Range is [0.011285, 0.460679]

FSD3

0.2

FSD3 Range is [0.002186, 0.337878]

After user clicks SEARCH, application will return results in popup window.

Result page will contain:

1. Show counts of matched records and slides.

Project 4

AboutServicesContact

Search Result

Statistic info

Find 24 records in 13 slides.

Please select a feature to view statistic info:

View statistic info

Please select slide to view:

TCGA-06-0216-01Z-00-DX3:004096.008192

View the slide

2. Select any feature to view histogram.

Project 4

AboutServicesContact

Search Result

Statistic info

Find 24 records in 13 slides.

Please select a feature to view statistic info:

FSD1

Area
Perimeter
Eccentricity
Circularity
MajorAxisLength
MinorAxisLength
Extent
Solidity
FSD1
FSD2
FSD3
FSD4
FSD5
FSD6
HematoxylinMeanIntensity
HematoxylinMeanMedianDifferenceIntensity
HematoxylinMaxIntensity
HematoxylinMinIntensity
HematoxylinStdIntensity

Please select slide to view:


TCGA-06-0216-01Z-00-DX3:000000.012288

View the slide

Please select a feature to view statistic info:

MajorAxisLength

View statistic info



3. Select any result slide. It will show preview of the image (200*200 pixels) along with detail information of the cells which match the query in the specific slide.

Project 4

About Services Contact

Please select slide to view:

TCGA-06-0216-01Z-00-DX3:032768.008192

View the slide

TCGA-06-0216-01Z-00-DX3:004096.008192
TCGA-06-0216-01Z-00-DX3:012288.028672
TCGA-06-0216-01Z-00-DX3:028672.024576
TCGA-06-0216-01Z-00-DX3:028672.008192
TCGA-06-0216-01Z-00-DX3:008192.004096
TCGA-06-0216-01Z-00-DX3:024576.012288
TCGA-06-0216-01Z-00-DX3:032768.024576
TCGA-06-0216-01Z-00-DX3:016384.008192
TCGA-06-0216-01Z-00-DX3:004096.012288
TCGA-06-0216-01Z-00-DX3:008192.008192
TCGA-06-0216-01Z-00-DX3:016384.028672
TCGA-06-0216-01Z-00-DX3:000000.012288
TCGA-06-0216-01Z-00-DX3:032768.008192

Cell in this slide:


Number	Details
1	Area=24.000000;Center=00;Circularity=0.843052; CytoplasmEnergy=0.008294; CytoplasmEnergyGradMag=0.020549; CytoplasmEntropy=7.087340; CytoplasmEntropyGradMag=5.893606; CytoplasmKurtosis=1.965312; CytoplasmKurtosisGradMag=6.217280; CytoplasmMaxIntensity=255.000000; CytoplasmMeanCanny=0.221289; CytoplasmMeanGradMag=28.995674; CytoplasmMeanIntensity=153.000810; CytoplasmMeanMedianDifferenceIntensity=8.374391; CytoplasmMinIntensity=45.930630; CytoplasmSkewness=0.233608; CytoplasmSkewnessGradMag=1.598040;

Please select slide to view:

TCGA-06-0216-01Z-00-DX3:016384.008192

[View the slide](#)

Slide image:



Cell in this slide:


Number	Details
1	Area=101.000000;Center=00;Circularity=0.502322; CytoplasmEnergy=0.013041; CytoplasmEnergyGradMag=0.026912; CytoplasmEntropy=6.667303; CytoplasmEntropyGradMag=5.614083; CytoplasmKurtosis=6.694879; CytoplasmKurtosisGradMag=6.951428; CytoplasmMaxIntensity=255.000000; CytoplasmMeanCanny=0.228319; CytoplasmMeanGradMag=21.428040; CytoplasmMeanIntensity=80.723768; CytoplasmMeanMedianDifferenceIntensity=12.079175; CytoplasmMinIntensity=15.891138; CytoplasmSkewness=1.887855; CytoplasmSkewnessGradMag=2.011627;

Please select slide to view:

TCGA-06-0216-01Z-00-DX3:004096.012288

[View the slide](#)

Slide image:



Cell in this slide:

Number	Details
1	Area=27.000000;Center=00;Circularity=0.618956; CytoplasmEnergy=0.016495; CytoplasmEnergyGradMag=0.030682; CytoplasmEntropy=6.139950; CytoplasmEntropyGradMag=5.260883; CytoplasmKurtosis=2.965653; CytoplasmKurtosisGradMag=4.376666; CytoplasmMaxIntensity=156.727509; CytoplasmMeanCanny=0.262009; CytoplasmMeanGradMag=21.109020; CytoplasmMeanIntensity=79.677201; CytoplasmMeanMedianDifferenceIntensity=1.078881; CytoplasmMinIntensity=22.936062; CytoplasmSkewness=0.239667; CytoplasmSkewnessGradMag=0.999400;

Project 4

About
Services
Contact

Please select slide to view:

TCGA-06-0216-01Z-00-DX3:000000.012288

View the slide

Slide image:

Cell in this slide:

Number	Details
1	Area=24.000000;Center=00;Circularity=0.843052;CytoplasmEnergy=0.008294;CytoplasmEnergyGradMag=0.020549;CytoplasmEntropy=7.087340;CytoplasmEntropyGradMag=5.893606;CytoplasmKurtosis=1.965312;CytoplasmKurtosisGradMag=6.217280;CytoplasmMaxIntensity=255.000000;CytoplasmMeanCanny=0.221289;CytoplasmMeanGradMag=28.995674;CytoplasmMeanIntensity=153.000810;CytoplasmMeanMedianDifferenceIntensity=8.374391;CytoplasmMinIntensity=45.930630;CytoplasmSkewness=0.233608;CytoplasmSkewnessGradMag=1.598040;
6	Area=115.000000;Center=00;Circularity=0.393893;CytoplasmEnergy=0.017078;CytoplasmEnergyGradMag=0.033073;CytoplasmEntropy=6.196777;CytoplasmEntropyGradMag=5.220700;CytoplasmKurtosis=6.278370;CytoplasmKurtosisGradMag=7.688766;CytoplasmMaxIntensity=184.696899;CytoplasmMeanCanny=0.240741;CytoplasmMeanGradMag=16.884085;CytoplasmMeanIntensity=62.908710;CytoplasmMeanMedianDifferenceIntensity=5.599668;CytoplasmMinIntensity=4.881145;CytoplasmSkewness=1.468701;CytoplasmSkewnessGradMag=1.789895;CytoplasmStdGradMag=12.224514;CytoplasmStdIntensity=24.680331;CytoplasmSumCanny=130.000000;Eccentricity=0.962249;Extent=0.357143;FSD1=0.099506;FSD2=0.117249;FSD3=0.159859;FSD4=0.225776;FSD5=0.055435;FSD6=0.038448;HematoxylinEnergy=0.147524;HematoxylinEnergyGradMag=0.022155;HematoxylinEntropy=4.577323;HematoxylinEntropyGradMag=5.805066;

- If user want to view original size slide with highlighted cells (4096*4096 pixels), Download is available. Click Download and Save File. Downloaded image is in the folder UserHome_Directory/Downloads/

